

A collection of TePe dental hygiene products is displayed on a teal background. On the left, three toothbrushes are arranged vertically: a blue one, a yellow one, and an orange one. In the center, there is a roll of white dental floss with a blue label that reads 'TePe DENTAL FLOSS 40 cm'. To the right of the floss, there are three more toothbrushes: a purple one labeled 'SELECT X-SOFT', a pink one with a unique handle design, and a light blue one labeled 'COMPACT X-SOFT'. At the bottom right, there is a small yellow toothbrush cleaner and a small orange object.



Table of contents

Tiny things for big smiles	3
Prevention	4
Sustainability	5
Interdental cleaning	6
Toothbrush & toothpaste	14
Dry mouth	24
Speciality brushes	28
Healthy smiles	32
Education	38
Customised products	41
Marketing material	42



Some products have symbols indicating their suggested use and available packaging types. Here is a reference guide for all symbols:

Areas of use:



Basic oral care for those with healthy mouths



Suitable for children and/or juniors



Suitable for those with braces



Suitable for those with implants



Suitable for those with gingivitis/periodontitis

Packaging:



Flow pack



Blister



Bag



Carton pack



Box



Tube



Bottle

Reference:

1. Greenhouse gas emissions are calculated in accordance with the GHG Protocol and ISO 14040 and ISO 14044.

2. TePe sources 100% raw materials for the toothbrush- and interdental brush handles, as well as the flosser, tongue cleaner, travel case, cap and extra grip, from either renewable raw materials or certified renewable raw materials via the mass balance approach, the latter independently certified by the

International Sustainability and Carbon Certificate PLUS (ISCC PLUS).

3. Furness S, et al. Interventions for the management of dry mouth: topical therapies. Cochrane Database Syst Rev. 2011 Dec 7;(12):CD008934.

4. Key findings within the group saying dry mouth is affecting their daily life (13%). TePe – Mouth Dryness in UK, web panel questionnaire study with 1000 respondents, Transvector, April 2023.

Tiny things for big smiles

It's often the tiny things in life that make it meaningful, like starting the day by sharing a smile with a friend, colleague or a stranger. The same is true for our oral health; a few minutes of toothbrushing, flossing, and in-between brushing can greatly affect our overall health and wellbeing.

It's a tiny effort, but a huge deal.

Taking care of your mouth means taking care of your health. That's why we continue innovating, educating, and providing a broad range of products for better oral health.

We started in 1965, and have continued to find clever solutions ever since. Look through this catalogue and discover the tiny things TePe offers to make people lead healthier lives with healthy smiles.



Prevention starts with daily care

You have probably heard us talk about prevention many times. That's because it's really important to us. Caring for your teeth and gums is more than just a routine to make your mouth feel fresh; it's essential for preventing common issues like cavities and gum disease. In other words, brushing, flossing, and regular

use of mouthwash – can save you from more extensive and expensive treatments in the future. Our products are designed to make daily prevention simple, effective, and enjoyable. Investing a little bit in our oral health today, can contribute to a healthier tomorrow.

1.



Education about prevention

TePe raises awareness about the positive impact of prevention with the aim to...

- Reduce development of oral diseases
- Reduce reoccurrence of oral diseases
- Improve or maintain level of oral diseases

2.

Increased knowledge improves health

Preventive care leads to healthier mouths which supports general health through...

- Functional wellbeing
- Psychological and social wellbeing
- The absence of discomfort or pain



3.



Improved health is good for the planet

Good general and oral health has a positive impact on the environment since it can lead to...

- Fewer resources used
- Less waste generated
- Reduced emissions

Sustainability is the healthy approach

Taking care of your teeth is great for your smile. Did you know that it can also contribute to creating a healthier planet? Healthy teeth often mean a lower environmental impact, as it can reduce the need for treatments and the emissions associated with dental repairs.

What about our products and production? In Sweden, all of our products are made using 100% renewable energy, boosted by our own solar park. We actively work to reduce greenhouse gas emissions and are phasing out virgin fossil plastics in favour of renewable options.^{1,2} Sustainability is a journey where improvements are made step by step, where every action counts – both big and small. That is what we call a healthy approach. Take a look at some of our milestones.

This is the beginning of our sustainability journey. TePe is founded with the launch of the triangular wooden dental stick, which is made of Swedish wood and developed in collaboration with the Faculty of Odontology in Malmö, Sweden.

1965




1993

TePe starts using hydro power. The production site is now powered by 100% renewable energy.

Innovation of preventive care: the TePe Interdental brush is launched.



2016




2019



A solar panel park is installed to power our production.

2018



TePe's cross-functional sustainability team is established.

Honeybees find a new home on our roof.

2022

All the handles for the TePe tooth-brushes and interdental brushes are from now on made with renewable or mass balance materials.



2021

Linking 7 of the Global Goals to our strategy.



2020

2030

Contribute to 25 million healthy smiles.



2023

Launching a tooth-brush with wooden handle and replaceable brush head and a dental floss made from recycled plastic bottles.



2024

Greenhouse gas mapping of scope 1, 2 and 3.



2025

100% of the packaging should now include renewable, recycled or FSC-certified materials.




Interdental cleaning

A special focus on interdental cleaning will keep those gaps gleaming and is essential for maintaining a healthy smile. We offer a range of interdental products designed to suit every individual and every type of gap.

Recommend the correct interdental cleaning product

When there are enough spaces between teeth, an interdental brush is the ideal tool. It can effectively reach and clean the nooks and crannies of the teeth, including the concave surfaces that floss may miss.

In narrow spaces where the gums completely fill the gap between the teeth, using floss or a toothpick is often sufficient to keep the area clean.






TePe® Interdental Brush Original

According to studies, interdental brushes are the most effective tool for removing plaque between teeth. Our brushes are developed together with dental experts, to ensure the highest functionality, comfort and quality. We offer 9 colour-coded ISO standard sizes, so there is a size for every type of gap, both big and small. Matching the gaps with the right colour makes it easy to remember what size you need. All the handles of our interdental brushes are made with renewable materials through the mass balance approach.²


The bristles cover the entire wire for safe and gentle cleaning.

Plastic-coated wire and rounded tip for safety and comfort.



The smaller brush sizes (0-3) have a flexible neck for increased durability and access.

The larger brush sizes have a thicker, bendable wire.



Ergonomic handle for a stable grip.



It is user-friendly thanks to the colour-coded sizes.



The plastic cap can be used for hygienic protection or to extend the handle.



TePe® Interdental Brush Original is also available in a mixed pack.




112 299




TePe® Interdental Brush Original


Colour: Pink
ISO size: 0
PHD size: 0.6 mm
Wire size: 0.4 mm




Colour: Orange
ISO size: 1
PHD size: 0.7 mm
Wire size: 0.45 mm




Colour: Red
ISO size: 2
PHD size: 0.9 mm
Wire size: 0.5 mm




Colour: Blue
ISO size: 3
PHD size: 1.1 mm
Wire size: 0.6 mm




Colour: Yellow
ISO size: 4
PHD size: 1.3 mm
Wire size: 0.7 mm




Colour: Green
ISO size: 5
PHD size: 1.6 mm
Wire size: 0.8 mm




Colour: Purple
ISO size: 6
PHD size: 2.0 mm
Wire size: 1.1 mm



Colour: Grey
ISO size: 7
PHD size: 2.6 mm
Wire size: 1.3 mm



Colour: Black
ISO size: 8
PHD size: 3.4 mm
Wire size: 1.5 mm




TePe® Interdental Brush Extra Soft


- Extra soft bristles for gentle and safe cleaning.
- Six sizes, all with plastic-coated wire.
- For patients with sensitive teeth, sore gingiva, dry mouth, or after oral surgery*.




Colour: Orange
ISO size: 1
PHD size: 0.7 mm
Wire size: 0.45 mm




Colour: Red
ISO size: 2
PHD size: 0.9 mm
Wire size: 0.5 mm




Colour: Blue
ISO size: 3
PHD size: 1.1 mm
Wire size: 0.6 mm




Colour: Yellow
ISO size: 4
PHD size: 1.3 mm
Wire size: 0.7 mm




Colour: Green
ISO size: 5
PHD size: 1.6 mm
Wire size: 0.8 mm



Colour: Purple
ISO size: 6
PHD size: 2.0 mm
Wire size: 1.1 mm



* A dental professional should always be consulted regarding brushing after oral surgery.



TePe Interdental Brush
Extra Soft is also available
in a mixed pack.


Colour: Pink
ISO size: 0
PHD size: 0.6 mm
Wire size: 0.4 mm



154 610 154 410



Colour: Orange
ISO size: 1
PHD size: 0.7 mm
Wire size: 0.45 mm



154 620 154 420



Colour: Red
ISO size: 2
PHD size: 0.9 mm
Wire size: 0.5 mm



154 630 154 430



Colour: Blue
ISO size: 3
PHD size: 1.1 mm
Wire size: 0.6 mm



154 640 154 440




Colour: Yellow
ISO size: 4
PHD size: 1.3 mm
Wire size: 0.7 mm



154 650 154 450



Colour: Green
ISO size: 5
PHD size: 1.6 mm
Wire size: 0.8 mm



154 660 154 460



- Plastic-coated wire for safe use.
- Long, flat handle for a stable grip.
- Angled brush head for easier access, especially between the posterior teeth.



154 699

TePe Angle™ is also available in a mixed pack.



TePe EasyPick™ – easy cleaning in between



This soft silicone toothpick efficiently removes plaque and pesky food debris to keep your teeth and gums healthy and clean. It's flexible and comfortable to use, even if you have braces, implants, or crowns.

Unique design and silicone coating ensure efficient and comfortable cleaning.

Rounded tip for comfortable use.

The unique lamellae provide a wide contact area with the surfaces between the teeth.


Overlapping conical sizes to fit narrow as well as wide spaces between teeth.



A handy pocket case for TePe EasyPick™ is included for easy cleaning on-the-go.

TePe EasyPick™

This product is a good choice for healthy mouths, in the case of gingivitis, and as an on-the-go alternative for everyone.



Colour: Orange

Size: Extra Small/Small



242 310



242 320

Colour: Blue

Size: Medium/Large



242 350



242 360

Colour: Purple


Size: Extra Large



242 330

TePe Mini Flosser™

- Pre-loaded dental floss holder for easy and convenient use in narrow spaces.
- Strong, slim and shred-resistant floss.
- Bite plane for safe and gentle insertion.
- Handle made with renewable materials through the mass balance approach.²



Insert the floss by gently biting onto the bite plane.



622 556



622 156

TePe® Gingival Gel

Size: 20 ml

- Contains chlorhexidine (0.2%) to fight bacteria, and fluoride (0.32% NaF) to strengthen the tooth enamel and protect sensitive root surfaces.
- Non-abrasive, without alcohol and colourants.
- Fresh mint flavour and smooth formula.
- Provides extra caries protection.



982 373



Use the cap as a bowl to apply gel.



TePe® Travel Case

- Practical travel case for TePe® Interdental Brush, TePe Mini Flosser™ or TePe EasyPick™.
- Air holes to allow bristles to dry.
- For hygienic storage during travel.
- Made with renewable materials through the mass balance approach.²



922 028



TePe® Dental Floss

Length: 5 m / 40 m


- Floss is recommended for tight interdental spaces, where an interdental brush will not fit.
- Expanding floss that suits different gaps as it adapts to the dental space.
- Slides easily between teeth due to its vegetable wax and avocado oil coating.
- Fresh mint flavour.



642 095
(5 m)



642 230
(40 m)




Gently move the floss along each tooth surface.



672 255



672 655



TePe® Bridge & Implant Floss

Amount: 30 pcs

- Wide and spongy floss for efficient cleaning.
- The short length and sturdy ends ease insertion and use.
- Invaluable for cleaning of implants, bridges and braces.
- Daily use contributes to healthy gums.





Made in
Sweden

TePe

**Daily™
Toothpaste**

Cavity prevention

Mild peppermint

With fluoride
Low-foaming
0% SLS
Vegan

Fluorid-Zahnpasta | Milder Minzgeschmack | Wenig Schaum
Fluoride tandpasta | Milde muntmaak | Laagschuimend
Dentifricio al fluoro | Aroma delicato di menta | Leggero
schiumoso
Dentifrice au Fluor | Saveur menthe douce | Peu moussant

TePe

TePe

Toothbrushes & toothpaste

High-quality, well-designed, effective and yet gentle on teeth and gums. TePe's toothbrushes come in various models, sizes, bristle types, and colours, making daily oral care easy and user-friendly for people of all ages.

And what would toothbrushes be without their toothpaste? We are proud to offer our own toothpaste range that perfectly complements our toothbrushes. The TePe Daily & Pure range offers mild fluoride toothpaste for every family member and people with sensitive or dry mouths.


A simple 3-step routine for a healthier smile

Always start by cleaning between your teeth to remove plaque. Using either an interdental brush or a dental floss is actually the most effective way to remove plaque between the teeth.



1.

The tongue cleaner is the magic wand against bad breath! The TePe Tongue Cleaner has three blades for a triple cleaning effect. It's easy to use, removes bacteria with just a few strokes and prevents bad breath.



2.

Start by brushing the insides of your teeth to make sure they get as much love as the outsides. Brush twice a day for two minutes for optimal plaque removal and fluoride distribution.




3.

TePe Supreme™ – supreme plaque removal with a gentle feel

Our innovative toothbrush that offers a longer reach, thanks to the two-level bristle design. The longer bristles clean further between the teeth and along the gum line, while the shorter bristles gently clean the rest of the tooth surfaces. Developed in collaboration with dental experts. And by the way, all the handles of our toothbrushes are made with renewable materials through the mass balance approach.²




Unique two-level
bristle design for
excellent reach.



Can be angled in hot
water for better access.

Easy-grip handle
with thumb pad for
extra comfort.

End-rounded and
tapered filaments.



TePe Supreme™

Size: Regular

Bristles: Soft

- Unique two-level bristle design for excellent reach.
- Easy-grip handle with thumb pad for extra comfort.



352 100



TePe Supreme™ Compact

Size: Compact

Bristles: Soft

- Unique two-level bristle design for excellent reach and easy-grip handle with thumb pad.
- Suitable for adolescents or for adults who prefer a smaller brush head.



392 120



TePe Nova™

Size: Regular

Bristles: Extra soft (pink), soft (yellow) & medium (blue)

- Long brush tip and end-rounded bristles for effective and gentle cleaning.
- Easy-grip handle with thumb pad for extra comfort.



312 110
Extra soft



312 120
Soft



312 130
Medium



The 4x2 rule is worth sharing



Use about 2 cm of fluoride toothpaste for adults.

(Children 0–2 years: a rice grain-sized amount.
3–6 years: a pea-sized amount)



Brush your teeth morning and evening.



Brush for at least two minutes.



Avoid eating or drinking for two hours to get the most benefit from the fluoride.

TePe Select™ Regular

Size: Regular

Bristles: Extra soft, soft & medium

- Tapered brush head for easy access.
- End-rounded bristles for gentle and effective cleaning.
- User-friendly handle design.



322 687
Extra soft



322 688
Soft



322 689
Medium



TePe Select™ Compact

Size: Compact

Bristles: Extra soft, soft & medium

- Small, tapered brush head for easy access.
- End-rounded bristles for gentle and effective cleaning.
- Suitable for children and adults who prefer a smaller brush head.



332 687
Extra soft



332 688
Soft



332 689
Medium



TePe Colour™

Size: Regular

Bristles: Soft

- Tapered brush head for easy access.
- End-rounded bristles for gentle and effective cleaning.
- User-friendly handle design.



322 683
(01-06)



322 383



TePe Colour™ Compact

Size: Compact

Bristles: Extra soft

- Small, tapered brush head for easy access.
- End-rounded bristles for gentle and effective cleaning.
- Suitable for children or for adults who prefer a smaller brush head.



332 682
(01-06)



332 482



TePe Mini™ Bolibompa

Size: Mini

Bristles: Extra soft

- Extra small, tapered brush head and end-rounded bristles for gentle and effective cleaning.
- Recommended from the first tooth.



301 686
Extra soft



301 386
Extra soft



11 186
tra soft



TePe Mini™

Size: Mini

Bristles: Extra soft & soft

- Extra small, tapered brush head and end-rounded bristles for gentle and effective cleaning.
- Recommended from the first tooth.
- Available with or without charming prints created by Swedish Isa Form.



382 687
Extra soft



382 310
Extra soft



New
design



TePe Kids™

Size: Compact

Bristles: Extra soft & soft

- Toothbrush with colourful print, suitable from the age of 6 years.
- Small, tapered brush head for easy access.
- End-rounded bristles for gentle and effective cleaning.



339 691
Extra soft



339 692
Soft



339 470
Extra soft



New
design



A good routine for healthy teeth in the future

Kids need help with brushing their teeth

Until they're 10-12 years or so.
Why not make it fun by adding a song or a story?



Create a habit

Brush the inside, outside and biting surfaces in the same order to make sure no side is forgotten. Remember the gum line where bacteria can easily build up and cause inflammation.





TePe Choice™ – brush away in a natural way

TePe Choice™ is an innovative and user-friendly toothbrush with a reusable handle made from Swedish beechwood, and three replaceable brush heads.

Size: Regular

Bristles: Soft

- Tapered brush head for easy access.
- Plant-based coating on the handle to minimise moisture absorption.
- Soft, end-rounded bristles for gentle and effective cleaning.
- Brush head made with renewable materials through the mass balance approach².



362 820



362 820-01
E-tail (green)




362 820-02
E-tail (pink)



TePe Classic™

Size: Regular
Bristles: Soft

- User-friendly handle design.
- End-rounded bristles for gentle and effective cleaning.
- For those who prefer a bigger brush head.



342 120
Soft




TePe PlaqSearch™

Size: 10 tablets/per pack


- A disclosing product for quick and simple self-check identification of plaque.
- Two-tone colour indicators.
- Teaching aid for better oral health.



992 680



992 677
(20 tablets)



Differentiates
new (red) and
old (blue) plaque.




TePe® Extra Grip

- Provides a comfortable, stable grip which helps patients with reduced manual dexterity.
- Fits the Select toothbrush and the specialty brushes.
- Easy to insert and remove by tapping the toothbrush handle on a hard surface.
- Weighs only 30 grams.
- Made with renewable materials through the mass balance approach.²



912 210



Insert/remove the
brush by tapping
the handle against
a table.



TePe® Brush Head Cap

- Small, practical travel cap.
- Air-holes ensure quick drying of the brush.
- Suits all TePe toothbrushes and special brushes (except the Denture and Classic specialty Brush).
- Made with renewable materials through the mass balance approach.²



922 024





TePe Daily™ Toothpaste

A high-quality toothpaste range that offers options for the whole family. It's effective against caries, gentle on the oral mucosa, and without unnecessary additives. TePe Daily is available in three varieties, for adults, kids, and babies, all with age-appropriate fluoride levels.

TePe Pure™ Toothpaste

The ideal choice for sensitive oral mucosa and dry mouth. The toothpaste come without foaming and colouring agents. Available in either a very mild taste or a flavourless option.



Good Environmental Choice

TePe has applied for the Swedish eco-label Bra Miljöval (Good Environmental Choice) and earned a certification for our toothpaste. We are proud to be the only toothpaste solid in Sweden with this environmental certification.

	TePe Daily™	TePe Pure™ (peppermint)	TePe Pure™ (unflavoured)
Caries prevention	✓✓✓	✓✓✓	✓✓✓
Sensitive mucosa		✓✓	✓✓✓
Dry mouth		✓✓✓	✓✓✓
Sensitive to SLS	✓✓✓	✓✓✓	✓✓✓
Sensory issues*		✓	✓✓✓
Flavour	Mild peppermint	Extra mild peppermint	Unflavoured
Foam	Low-foaming	No foam	No foam
Colour	White	No colour	No colour

* I.e. in patient care or people with developmental disabilities. Also suitable for children 7 years or above.

TePe Daily™ Baby

Sodium fluoride: 1000 ppm F⁻

Size: 50 ml, 15 ml

Age: 0-6 years

- Easier brushing without foam.
- No added flavouring or colouring agents.
- Vegan, SLS-free.



Use a rice grain sized amount of fluoride toothpaste for ages 0-2 years.



50 ml
802 112



15 ml
802 212



TePe Daily™

Sodium fluoride: 1450 ppm F⁻

Size: 75 ml, 15 ml

Age: From 7 years of age

- Low-foaming for comfortable brushing.
- Mild peppermint flavour.
- Vegan, SLS-free.



75 ml
802 110



15 ml
802 210



Use 2 cm of
fluoride tooth-pas-
te.



TePe Pure™ Mild peppermint

Sodium fluoride: 1450 ppm F⁻

Size: 75 ml, 15 ml

Age: From 7 years of age

- No foam for comfortable brushing.
- Extra mild peppermint flavour.
- No colouring agents.
- Vegan, SLS-free.



75 ml
812 114



15 ml
812 214



Use 2 cm of
fluoride tooth-pas-
te.



TePe Daily™ Kids

Sodium fluoride: 1000 ppm F⁻

Size: 75 ml, 15 ml

Age: 0-6 years

- Easier brushing without foam.
- Extra mild peppermint flavour.
- Vegan, SLS-free.



75 ml
802 111



15 ml
802 211



Use a pea sized
amount of fluoride
toothpaste for ages
3-6 years.



TePe Pure™ Unflavoured

Sodium fluoride: 1450 ppm F⁻

Size: 75 ml, 15 ml

Age: From 7 years of age

- No foam, for comfortable brushing.
- Unflavoured.
- No colouring agents.
- Vegan, SLS-free.



75 ml
812 113



15 ml
812 213



Use 2 cm of
fluoride tooth-pas-
te.



Dry mouth

Dry mouth is an often overlooked issue that affects 1 in 5 people.³ It is closely linked to oral health. We are proud to offer this important piece of the puzzle of improving oral health, with a wide range of products that soothe, provide comfort and protect teeth.

What do dry mouth sufferers say?⁴



26%

say dry mouth makes it
difficult to speak



33%

of those affected
in their daily life say
that their **self confidence**
is **negatively affected**

From those suffering
more severely,

35%

say it **affects**
their sleep



62%

do not talk about their
dry mouth condition

Certain
illnesses,
medications
and menopause
can cause
dry mouth



Now offering comfort for dry mouth

TePe's hydrate range is designed to provide easy and long-lasting comfort, wherever you are and whenever you need it. The products moisten, soothe and refresh.

The hydrating mouthwashes with fluoride also give additional caries protection. Due to the high fluoride in the mouthwash, we recommend using it 1-2 times/day, preferably not in connection to toothbrushing, for extra caries protection.

Unflavoured or flavoured?

We recommend unflavoured products for people with little or no saliva production and very sensitive or damaged mucosa.



TePe® Hydrating Mouthwash Mild apple peppermint

Sodium fluoride: 0,2%

Saliva-stimulating & moisturising: Malic Acid, Betaine and glycerin

Size: 500 ml

Age: From 12 years of age

- Moistens the mucosa and leaves a pleasant feeling in the mouth.
- Saliva-stimulating.
- Mild taste of apple and peppermint.
- Can be used twice a day.
- Practical dosing cup.



832 314



TePe® Hydrating Mouthwash Unflavoured

Sodium fluoride: 0,2%

Moisturising: Glycerin

Size: 500 ml

Age: From 12 years of age

- Moistens the mucosa and leaves a pleasant feeling in the mouth.
- Gentle and unflavoured.
- Contains no alcohol.
- For very sensitive mouths.
- Can be used twice a day.



832 313



TePe® Hydrating Mouth Gel Mild peppermint

Caries prevention: Xylitol

Moisturising: Betaine and glycerin

Size: 50 ml

Age: All ages

- Moistens, soothes and provides an immediate feeling of comfort.*
- Convenient, hygienic and easy to use also on-the-go.
- Can be used as often as needed, day and night
- Fresh feeling of mild peppermint.



842 314

*The duration of the comfort is individual; people with no or very little saliva usually experience a more prolonged effect.



TePe® Hydrating Mouth Gel Unflavoured

Xylitol

Moisturising: Betaine and glycerin

Size: 50 ml

Age: All ages

- Moistens, soothes and provides an immediate feeling of comfort.*
- Convenient, hygienic and easy to use also on-the-go.
- Can be used as often as needed, day and night.
- Unflavoured and gentle.



842 313



TePe® Hydrating Mouth Spray Mild apple peppermint

Caries prevention: SMFP 40 ppm F, Xylitol

Saliva-stimulating & moisturising: Malic Acid and glycerin

Size: 50 ml

Age: From 12 years of age

- Saliva-stimulating and gives an instant feeling of freshness.
- Convenient on-the-go.
- Can be used max 20 times per day. Suitable during the day.
- Mild apple and peppermint.



862 314



TePe® Lubricating Mouth Spray** Unflavoured

Moisturising: Sunflower oil

Size: 50 ml

Age: All ages

- A lubricating mouth spray for very dry mouths.
- Mild taste.
- Convenient on-the-go.
- Use as often as needed.



862 313

**Currently available in selected markets.



New

Speciality brushes

Wisdom teeth, dental implants, crowns/bridges or braces – no matter your needs, TePe can help make cleaning easier. TePe offers a wide range of high-quality cleaning devices with excellent functionality for those with specific oral hygiene needs.




TePe Special Care™

Size: Regular (red) & Compact (blue)

Bristles: Ultra soft

- For patients with extremely sore and delicate tissues, e.g. as a result of trauma, radiotherapy treatment or after oral surgery.
- 12,000 (Regular) and 9,000 (Compact) ultra soft bristles for delicate cleaning.
- First brush stage after oral surgery.*



* A dental professional should always be consulted regarding brushing after oral surgery.



TePe Gentle Care™

Size: Regular

Bristles: Super soft

- For patients with dry mouth, sensitive root surfaces or sensitive oral tissues.
- Super soft filaments for gentle cleaning.
- Second brush stage after oral surgery.*



462 150



TePe® Denture Care

Size: Regular

- Effective and safe cleaning of removable dentures.
- Brush head can be bent backwards for improved grip and accessibility.
- Designed to resemble a regular toothbrush.



422 368



Long, strong filaments for efficient cleaning of dentures.



TePe® Tongue Cleaner

- Three blades for a triple cleaning effect.
- Helps prevent bad breath (halitosis).
- Removes bacteria in a few strokes.
- Made with renewable materials through the mass balance approach.²



942 385



TePe Compact Tuft™

TePe Compact Tuft is an angled brush designed for precision. It features a small, dome-shaped brush with tightly packed short filaments. With its stable handle and dense, firm tuft, this brush makes it easier to clean even the hardest-to-reach areas. And by the way, all the handles of our specialty brushes are made with renewable materials through the mass balance approach.²

Dome-shaped brush head with dense, rounded tuft.

Tightly packed short bristles that give a firm feel.



TePe Compact Tuft™

- Optimal for precision cleaning.
- Dome-shaped brush head with dense, rounded tuft.
- Ideal for cleaning along the gum line, on the occlusal surfaces of emerging teeth, or impacted wisdom teeth.



432 100



TePe Universal Care™

- A unique angle provides easy access to the lingual and palatal surfaces of teeth and implants.
- Makes it easier to clean narrow and hard-to-reach areas.
- End-rounded bristles for gentle cleaning.



472 368



The angled neck makes it easier to clean the lingual surfaces.



TePe Interspace™

Bristles: Soft (yellow) & medium (blue)


- Ideal for cleaning furcations, implants, fixed braces or behind the last tooth.
- Brush tip can be attached to either side of the handle.
- Comes with 12 replaceable brush tips.



x12
412 364
Soft



x12
412 366
Medium



TePe® Implant Care Kit*

A collection of carefully selected maintenance products designed for effective plaque removal in and around implants, making the cleaning simple.

Products included:

- TePe Compact Tuft™
- TePe Implant Orthodontic™ Brush
- TePe Universal Care™
- TePe® Bridge & Implant Floss (30 pcs)
- TePe® Interdental Brushes Extra Soft mixed pack
- Patient Information Leaflet



962 687



TePe® Implant Orthodontic Brush

- Extra slim brush head for improved access.
- End-rounded bristles for gentle and effective cleaning.
- Optimal for cleaning buccal implant surfaces or orthodontic appliances.
- Long and thin neck for better access.



x1
442 100



Effective cleaning
around brackets
and wire.



TePe® Orthodontic Kit*

The kit features carefully selected self-care products and readily accessible tips and advice for orthodontic patients.

Products included:

- TePe Supreme™ Compact
- TePe Compact Tuft™
- TePe® Interdental Brush, ISO size 4 and 5
- TePe Angle™ Interdental brush, ISO size 4
- TePe EasyPick™, XS/S and M/L
- TePe® Orthodontic wax
- Patient Information Leaflet



962 686

*Currently available in selected markets.





Healthy smiles for all mouths

TePe offers a wide selection of products suitable for all stages of life – from the first teeth to special needs like oral care for braces and implants. You will find our product recommendations throughout the following pages.





Paedodontic care

To keep the smile bright and healthy for a lifetime, we suggest brushing as soon as that first tooth pops up. It is often recommended to start with interdental cleaning during the teenage years.



TePe Supreme™
Compact



TePe Compact
Tuft™



TePe Select™
Compact



TePe Mini™



TePe PlaqSearch™

The foundation of maintaining healthy gums is daily interdental cleaning and toothbrushing. Professionals play a key role in motivating patients and recommending the right products.


Daily care





Implant care

Getting the right advice on suitable products and cleaning techniques can make a big difference in keeping your implants healthy. Daily cleaning helps to prevent peri-implant diseases.



TePe® Interdental
Brushes



TePe® Implant
Orthodontic



TePe® Bridge
& Implant Floss




TePe Compact
Tuft™



TePe Universal
Care™



TePe®
Denture Care



TePe EasyPick™

Fixed orthodontic appliances, like braces, can quickly accumulate plaque around brackets, bands, and wires. To keep the smile healthy and bright, it's important to establish a good self-care routine that helps prevent gingivitis and caries..

Orthodontic care



TePe Universal
Care™




TePe Supreme™
Compact




TePe Interspace™



TePe PlaqSearch™




TePe® Interdental
Brushes



TePe® Implant
Orthodontic



TePe Compact
Tuft™



TePe EasyPick™



Want a bigger bite of science?

Sign up for TePe Clinical Update



Articles



Lectures



Webinars



Newsletter



Events



Sharing
materials



Four times a year, we will provide dental professionals and students with a bundle of inspiration to apply in clinical practice. We offer our newsletter Simply Science, a Webinar and Tips and Tricks.

Simply Science is a valuable resource keeping you updated with the latest scientific research in oral health. It is a brief summary of a recent scientific article, along with key insights that you can directly apply to your clinical

practice. Our Clinical and Education Team will focus on sharing topics related to TePe's core mission and the work of dental professionals - preventing oral disease and promoting healthy smiles worldwide.

We will continue with our highly appreciated webinars. Rest assured, as long as you register, you can access the recording even if you don't attend the live session. We don't want anyone to miss out on valuable knowledge.

Explore our knowledge
platform at tepe.com



Scan the
QR code
to sign up

Our global Clinical and Education Team:

RDH Lina Gassner Kanters

Clinical and Education Advisor

Dr Michaela von Geijer

Clinical and Education Advisor

Dr Sanjay Haryana

Clinical and Education Advisor





We provide clinicians with knowledge, tips, and tricks to achieve better compliance and treatment results.





Practice Box

- Practical storage of TePe's interdental brushes in the dental practice.
- 8 removable trays for IDB Original and Extra Soft, and 4 trays for TePe Angle. Each tray holds 15 brushes.
- The box includes a helpful chart that displays various brush sizes.

IDB Original /Extra Soft	TePe Angle™	Complete IDB, incl. 40 mix	Complete Angle, incl. 30 mix
			
922 045	922 041	2400022	2400023

TePe® Dental Model

- For patient demonstration and instruction.
- Shows crowding, various interdental spaces, furcation, an erupting wisdom tooth and tooth loss.
- Can be marked with a non-permanent pen for instructional use.
- Suitable for the product cases featured below.



952 051

Delivered in a brown hard paper box.



TePe® Product Demonstration Case

Prepacked with the majority of our products to show the variation of the product range. Includes extra space for marketing material and samples. Closes with a zipper. Carry it with ease via the shoulder strap or handle. Suited for sales representatives.



952 056



TePe® Product Presentation Case

Case featuring the major part of the TePe range. Suitable for dental clinics, for example.



952 054



Customised products

Put the clinic's logo and a tagline on a toothbrush and spread awareness about the services, while promoting better oral health. TePe's imprinting service is an excellent marketing opportunity for dental practices.



Logotypes, images and messages can be printed on our Select, Select Compact, Colour, Colour Compact and Mini toothbrushes. Choose between a variety of different handle colours and print colours for amazing designs.



Marketing & display material

TePe offers a range of information and demonstration materials to aid oral hygiene advice. Our product racks are perfect for display and storage.

Before ordering, please reach out to your TePe contact person, to make sure the design and product is available on your market.





Paper Bag

Item no: 256 INT



Patient leaflet Display

Item no: 931 050
Display with 3 pockets for patient folders.
Clear acrylic material.



Mini Display

Item no: 931 031
For 6 toothbrushes.
9 interdental brushes.
Clear acrylic material.



Toothbrush Display

Item no: 931 029
For 13 toothbrushes
(Select range).
Clear acrylic material.



Interdental Brush Display

Item no: 931 035
For 15 interdental brushes.
Clear acrylic material.



Floor Stand, Metal

Item no: 931 054
Floor stand on wheels.
21 hangers, 1 shelf.
2 leaflet holders.



Mini Flosser/Dental floss Dispenser Box

Item no: 622 570 (Mini flosser, 120 pcs)
642 571 (Dental floss, 80 pcs)
Dispenser box for Mini flosser /
Dental floss samples.



EasyPick Dispenser Box

Item no: 242 261 (XS/S)
242 262 (M/L)
242 263 (XL)
242 263 (Mixed XS/S, M/L)
Dispenser box for EasyPick
samples. 200 x 2-pack included.



EasyPick Display

Item no: 221INT
Cardboard table display.
For 3 SRPs, EasyPick.
(shelf ready packaging)



Sample Box

Item no: 152
Filled with a variety of product
samples and information.



Floor Stand, Carton

Item no: 932 092
Floor stand.
4 shelf.



Medium Table Display

Item no: 932 033
Table display in cardboard.
8 hangers.



Wall Display, Metal

Item no: 931 053
Wall display
12 hangers. 1 shelf.
(Plastic containers available
separately, item no. 931 038).



Table Display, Metal

Item no: 931 055
Table display with
rotating function.
8 hangers.



Table Display, Carton

Item no: 932 091
Fits 2 SRPs.
(Toothpaste holder available
separately, item no. 932 094).

*Please note the visual design might change.



TePe is certified according to the international quality standard ISO 9001 and the environmental standard ISO 14001.

TePe Oral Hygiene Products Ltd.

Eklund House | Cathedral Avenue

Wells | Somerset BA5 1FD

Tel: 01749 60 88 00 | Fax: 01749 60 87 79

info@tepe.co.uk | www.tepe.com

TePe Munhygienprodukter AB

Bronsåldersgatan 5 | 213 76 Malmö, Sweden

Tel +46 (0)40 670 11 00 | Fax +46 (0)40 670 11 11

info@tepe.com | www.tepe.com

All trademarks herein are the property of TePe Munhygienprodukter AB.