



TePe Supreme™/ Supreme™ Compact

Product Features

- Unique two-level filament design for excellent reach between the teeth and along the gum line
- Longer bristles and tapered brush-head for improved access
- Ergonomically designed non-slip handle for secure grip
- Neck can be angled without heating
- Available in regular or compact size
- For single person use only

Product Materials & Specifications

Handle Material: Polypropylene (ISCC Certified)

Grip Material: Thermoplastic elastomer (TPE).

Filament Material:

TePe Supreme = Polybutylene terephthalate (PBT) (~ 2400 filaments)

Tepe Supreme Compact = Polybutylene terephthalate (PBT) filaments (~1850 filaments)

Level of Softness: Soft

Colours: Pink, Orange, Yellow, Green Blue, Purple



Packaging Materials

Flowpacks: Polypropylene (PP).

Blister Pack-Blister: Polyethylene terephthalate (PET).

Blister Pack-Backcard: FSC cardboard.

Packaging Specifications

Flowpacks:

1pc/flowpack, 80 flowpacks/ white carton. Mixed colours.

Blister Packs:

1 pc/blister pack, 14 blister packs/ SRP carton. Mixed colours

3pcs/blister, 7 blisterpacks/ SRP carton. Mixed Colours

Recycling

Toothbrush: Household Waste

Flowpack: Container for Plastic

Blister Pack-Blister: Container for Plastic

Blister Pack-Backcard: Container for Corrugated Cardboard

Country of Origin

Sweden